

Product Line

Part Number	Description	Kit Price
Z-100	Feeding Bag for Liquids 5 - Individual bags each with feed tube and clamp all packaged in a zipper bag 1 - Product insert with instructions Use - Thin liquids, strained soups, supplements, medications and rinses	\$13.49
Z-200	Feeding Bag for Pureed Food 5 - Individual bags each with feed tube and clamp all packaged in a zipper bag 1 - Product insert with instructions Use - Thicker liquids, blended foods, thicker supplements (tubing is larger than Z-100)	\$14.95
Z-500	"Dinner Through a Straw" Recovery Book 1 - "Dinner Through a Straw" Recovery Book 1 - Product insert with instructions Use - Up-to-date information on the "how-to" of blending and maintaining adequate nutrition while on a liquid diet. Complete with recovery information and lots of great tasting "easy-to-make" recipes for liquid and blended meals.	\$7.25
Z-600	Patient Discharge Kit 3 - Z-100 Liquid Feeding Bag Kits 3 - Z-200 Pureed Food Bag Kits 1 - "Dinner Through a Straw" Recovery Book 1 - Product insert with instructions	\$83.25

Instructions for Use

1. Use Z-100 for thin liquids and the Z-200 for thicker products. The valve and feed tube are larger on the Z-200 in order to permit easier flow.
2. Attach feed tube to valve stem by pushing and turning feed tube until valve stem is fully inserted into tube. Attach clamp onto feed tube and pinch closed.
3. Open the resealable zipper top to fill with fluid. Fill to indicated line only. Use a folded paper plate or funnel when filling with pureed food. Pinch zipper close on one end and slide thumb and forefinger along zipper to seal. Repeat closing process.
4. Insert feed tube into mouth, release clamp and gently squeeze to dispense product.
5. Use warm or cold fluids only. Do not fill with hot liquids. Do not lift or carry bags by the feed tube.



800-824-0817

www.com-pac.com
www.bagstoreusa.com



For Easy Oral Care and Feeding



Resealable Bags
for
Liquid Nutrition & Oral Hygiene

"Healthier, happier patients!"





For Easy Oral Care and Feeding



Zip-N-Squeeze bags were designed for convenience and comfort for patients with facial and mouth surgeries who face difficulties with, or are unable to use, cups and straws for feeding. Designed and manufactured specifically to be

squeezed with no sucking required prevents pain, pressure and new incisions from bleeding. Gives jaws a rest and diverts liquid foods away from painful areas.

There is a **Zip-N-Squeeze** bag for liquid hydration, nutrition, medication and mouth care as well as one for pureed foods. **Zip-N-Squeeze** bags feature a soft, flexible tube which allows patients to gently irrigate difficult to reach areas where a toothbrush temporarily cannot be used. This provides a convenient method for frequent mouth care - you can mix whatever solution you prefer to use.

Swelling and pain can make it difficult for patients to drink from a cup, and discharge instructions usually state that patients should avoid using straws. **Zip-N-Squeeze** bags work wonderfully when filled with broth, liquids, protein drinks, medication, mouth wash and pureed foods. The pleated bottom allows the bag to stand upright, the pinch clamp prevents spills and the patented leakproof zipper will not rupture when pressure is applied.

Benefits:

- ✓ *Promotes hydration and nutrition intake for swollen and sore mouths and faces*
- ✓ *Easy to use, no mess*
- ✓ *Gentle squeeze of the bags dispenses contents - "No Sucking Required"*
- ✓ *Promotes better nutrition when cups and straws are difficult to use*
- ✓ *Leakproof resealable zipper stays closed when squeezing*
- ✓ *Soft, flexible tubing promotes increased comfort*
- ✓ *Lowers anxiety at meal time thereby increasing appetite*
- ✓ *Food can be pre-measured, filled in advance and stored in refrigerator or freezer*
- ✓ *Disposable - Reduces risk of bacterial contamination*
- ✓ *All materials are FDA approved for direct contact with food and liquids*

Zip-N-Squeeze products come highly recommended by the patients who use them.

Here's what they told us:

"They are so easy to use, I don't know how I would have gotten along without them."

"I was able to get more fluids in faster and be satisfied quicker."

"The syringe was so difficult to use, especially because I didn't feel well. When I got the Zip-N-Squeeze, it was the answer to my prayers, it's so easy and works so well."

"The Zip-N-Squeeze was so easy and held so much, the syringe was too much work and only held two tablespoons."

"Everyone should be able to use these bags. They made such a difference in my recovery!"

"Thank you for inventing such a wonderful and useful product."

And here's what their doctors told us:

"All our patients are to receive Zip-N-Squeeze kits before discharge. We've found them to make a real difference in their recovery. Patients actually take in more fluids because it's so easy to use."

"These bags definitely make hydration easier."

"All orthognathic patients would benefit from this product."